

Instructions for Use

ICELOCK® 125 RATCHET

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
EL	Οδηγίες Χρήσης	37
FI	Käyttöohjeet	41
NL	Gebruiksaanwijzing	45
PT	Instruções de Utilização	49
PL	Instrukcja użytkowania	53
TR	Kullanım Talimatları	57
RU	Инструкция по использованию	61
JA	取扱説明書	65
ZH	中文说明书	69
KO	사용 설명서	73

1

ENGLISH

Medical Device

DESCRIPTION

The device is a mechanical lock used to connect and release a locking liner into a prosthetic socket.

Components (Fig. 1)

1. Lock Body
2. Puck
3. Screws
4. Pin
5. Attachment Pin

INTENDED USE

The device is intended to connect and release to a prosthetic system that replaces a missing limb

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.

The weight limit for the device is 120 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur components and replacement parts in combination with this device. Use of other components may affect safety and performance of the device.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Align the Lock Body in the Puck. Apply medium strength threadlocker on the four Screws and torque them to 10 Nm.

Release unit

Shorten the release tube and Pin equally to the required length. Insert the Pin into the tube and screw the tube firmly and completely into the hole.

Attachment Pin

Apply Loctite 410/411 onto the threads on the Attachment Pin and screw it into the Locking Liner. Torque the Attachment Pin to 4 Nm.

Warning: Ensure that the pin is fully engaged. Extended use of the device without full engagement may affect safety and performance of the device.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

Warning: Inspect the lock and pin for wear at least once per year or if there is a change or loss in functionality. In case of wear, replace both lock and pin. Excessive wear may cause loss of suspension and an increased risk of falling.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P5 - 120 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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