
Instructions for Use

ICELOCK® 125 RATCHET
FABRICATION KIT WITHOUT
BRAID

EN		Instructions for Use	4
DE		Gebrauchsanweisung	8
FR		Notice d'utilisation	13
ES		Instrucciones para el uso	18
IT		Istruzioni per l'uso	23
NO		Bruksanvisning	28
DA		Brugsanvisning	32
SV		Bruksanvisning	36
EL		Οδηγίες Χρήσης	40
FI		Käyttöohjeet	45
NL		Gebruiksaanwijzing	49
PT		Instruções de Utilização	54
PL		Instrukcja użytkowania	59
CS		Návod k použití	64
TR		Kullanım Talimatları	68
RU		Инструкция по использованию	72
JA		取扱説明書	77
ZH		中文说明书	81
KO		사용 설명서	85

1

ENGLISH

Medical Device

DESCRIPTION

The device is a mechanical lock used to connect and release a locking liner into a prosthetic socket.

Components (Fig. 1)

1. Distal Spacer Pad
2. PVA Bag
3. Puck
4. Washer
5. Smart Cast Tool
6. Bolt
7. Tie Off
8. Proximal Groove

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.
The weight limit for the device is 120 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

1. Trim the Distal Spacer Pad to the size of the distal end of the cast, put it in place, and pull the inner PVA Bag over.
Note: If using Pressure Pads, adhere them onto the PVA Bag over the bony areas. Make sure that the PVA bag is completely dry.
2. Fill the heads of the four screws in the Puck with clay.
3. Position the Puck on top of the cast pulling the PVA Bag through the central hole.
4. Place the Washer from the SmartCast Tool inside the PVA Bag. Secure the Puck to the SmartCast Tool with the Bolt. Make sure that the Puck is orientated with the screws in AP/ML position.
5. Tie Off the PVA Bag and evacuate the air.
6. Fill the Puck with clay and seal it with sealing resin.
7. Tie on the lay-up material. Tie the first layer in the most Proximal Groove and the subsequent layers gradually more distally.
8. Apply the outer PVA Bag and place the aluminum ring inside it. Keep the ring 1–2 cm above the Puck.

9. Pour resin, allowing it to impregnate around the Puck.
10. Push the aluminum ring down over the Puck.
11. Finish lamination.
12. Remove the bolt and the socket from the cast. Trim and smooth edges.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning. Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY


Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3–5 years of use.

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Össur UK Ltd