


Instructions for Use

ICELOCK® 211 CLUTCH

| | | | |
|----|--|-----------------------------|-----|
| EN | | Instructions for Use | 5 |
| DE | | Gebrauchsanweisung | 10 |
| FR | | Notice d'utilisation | 16 |
| ES | | Instrucciones para el uso | 22 |
| IT | | Istruzioni per l'uso | 28 |
| NO | | Bruksanvisning | 34 |
| DA | | Brugsanvisning | 39 |
| SV | | Bruksanvisning | 44 |
| EL | | Οδηγίες Χρήσης | 49 |
| FI | | Käyttöohjeet | 55 |
| NL | | Gebruiksaanwijzing | 60 |
| PT | | Instruções de Utilização | 66 |
| PL | | Instrukcja użytkowania | 72 |
| CS | | Návod k použití | 78 |
| TR | | Kullanım Talimatları | 83 |
| RU | | Инструкция по использованию | 88 |
| JA | | 取扱説明書 | 94 |
| ZH | | 中文说明书 | 99 |
| KO | | 사용 설명서 | 104 |

1

ENGLISH

Medical Device

DESCRIPTION

The device is a mechanical lock used to connect and release a locking liner into a prosthetic socket.

Components (Fig. 1)

1. Assembly Bolt
2. Pin Tunnel
3. Lock Body
4. Dummy Screw
5. Silicone Seal
6. Delrin Guide
7. Clutch Mechanism
8. Push Button
9. Plaster Cast
10. T-Wrench
11. Attachment Pin

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Lamination Instructions

1. Flatten the distal end of Plaster Cast to center Delrin Guide. Attach with nails provided or double-sided tape maintaining original pin alignment.
2. Apply inner PVA bag and tie it around the neck (G) on Delrin Guide.
3. Trim away excess PVA and centrally locate Silicone Seal.
4. Place Lock Body onto Delrin Guide, ensuring its rotation allows for correct positioning of Clutch Mechanism.
5. Secure Lock Body and Pin Tunnel with Assembly Bolt (making sure that o-ring is between bolt / Pin Tunnel and Pin Tunnel / Lock Body).
Note: Be sure to tighten so that resin will not leak in between components.
6. Place some clay in the hole on bolt (A) and groove (B) on Dummy Screw.

7. Tie half of the lamination material into the groove on Lock Body (C) and remainder over the end of Lock Body (D) leaving Pin Tunnel and bolt head (E) exposed.
8. Apply outer PVA bag and finish lamination.
9. When curing has finished, remove Bolt, expose but do not remove Dummy Screw.
10. Before the second lamination ensure that Pin Tunnel is sealed.
11. After completion of the second lamination, expose and remove Dummy Screw.
12. Insert and tighten Clutch Mechanism to 15 Nm.

Warning: Make sure to tighten the Clutch Mechanism to 15 Nm.

Note: Push Button may be shortened by unscrewing from Lock Body and cutting it with hacksaw.

Thermoplastic Instructions

1. Flatten the distal end of Plaster Cast to center Delrin Guide. Attach with nails provided or double-sided tape maintaining original pin alignment.
2. Place Lock Body onto Delrin Guide, ensuring its rotation allows for correct positioning of Clutch Mechanism.
3. Secure Lock Body and Pin Tunnel with Assembly Bolt.
4. Vacuum form thermoplastic, ensuring that plastic conforms well to the Lock Body and Pin Tunnel.
5. Once the plastic has cooled, trim around Pin Tunnel and Dummy Screw.
Note: Before making the outer socket, always ensure that Assembly Bolt has been removed and Pin Tunnel resealed.
6. After manufacturing of outer socket, expose and remove Dummy Screw.
7. Insert and tighten the Clutch Mechanism to 15 Nm.

Warning: Make sure to tighten the Clutch Mechanism to 15 Nm.

Note: Push Button may be shortened by unscrewing from Lock Body and cutting it with hacksaw.

Attachment Pin

Apply Loctite 410/411 onto the threads on the Attachment Pin and screw it into the Locking Liner. Torque the Attachment Pin to 4 Nm.

Caution: Use Attachment Pins for Icelock 200 series only.

T-wrench supplied is specially designed to wind Push Button on Clutch Mechanism when donning socket.

Note: Use of any other tool may result in over-tightening.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning. Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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Össur UK Ltd