













Instructions for Use





ICELOCK® 214 CLUTCH AK OFFSET

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	12
FR	Notice d'utilisation	20
ES	Instrucciones para el uso	27
IT	Istruzioni per l'uso	34
NO	Bruksanvisning	41
DA	Brugsanvisning	47
SV	Bruksanvisning	54
EL	Οδηγίες Χρήσης	60
FI	Käyttöohjeet	67
NL	Gebruiksaanwijzing	74
PT	Instruções de Utilização	81
PL	Instrukcja użytkowania	88
TR	Kullanım Talimatları	95
RU	Инструкция по использованию	102
JA	取扱説明書	110
ZH	中文说明书	116
KO	사용 설명서	122

ENGLISH

Medical Device

DESCRIPTION

The device is a mechanical lock used to connect and release a locking liner into a prosthetic socket.

Components:

- A. Manufacturing Dummy
- B. Four Screws
- C. Silicone Ring
- D. Lock Body
- E. 4-Hole Gasket
- F. 4-Hole Dummy
- G. 4-Hole Screw
- H. Black Dummy Screw
- I. Grey Dummy Screw
- J. Central Screw
- K. Plastic Mesh
- L. Clutch Pin Cover
- M. Clutch Key
- N. Clutch Mechanism
- O. Clutch Pin

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.
The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.

The weight limit for the device is 136kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Prepare plaster negative as usual.

Preparing the Plaster Positive

1. Model the plaster positive to agree with the residual limb. Mold the distal end of the liner as precisely as possible.
2. Flatten the distal end with a plaster rasp until the diameter is exactly the same size as the Manufacturing Dummy (**Fig. 1a**).
3. For thermoplastic sockets only: Use a knife to carve out air channels on the flattened distal end to allow for airflow (**Fig. 2a**).
4. Use one of the two following methods to ensure airflow to the distal connector at the distal end:

- Drill a hole through the whole model from the distal center to the proximal end before fixing the Manufacturing Dummy (**Fig. 1a**).

OR

- Apply a thin stocking to the plaster.
5. Use the Four Screws to attach the Manufacturing Dummy to the plaster model (**Fig. 3a**).

Note: If using a moist plaster model, isolate properly.

Thermoplastic Sockets

1. Use the Central Screw to secure the Lock Body with Silicone Ring on the Manufacturing Dummy and the plaster model. Use the 4-Hole Screws to secure the 4-Hole Dummy (**Fig. 4a**).
2. Be careful with the orientation of the Clutch Mechanism. See **figure 5a** for the correct positioning of the black and grey dummies.
Note: The dummies have opposite threads.
3. Select the applicable material depending on the expected stresses and drape as usual (**Fig. 6a**).
4. Tightly wind a strap into the groove of the Lock Body during draping.
5. After fully hardening, grind and remove the material at the distal end until completely exposing the head of the Central Screw (**Fig. 7**).
6. Remove the Central Screw.
7. Grind carefully at the location of the 4-Hole Dummy until completely exposing the head of the 4-Hole Screws and 4-Hole Dummy (**Fig. 8**).
8. Remove the 4-Hole Screws and 4-Hole Dummy.
9. Grind carefully at the future location of the Clutch Mechanism until completely exposing the head of the Dummy Screw (**Fig. 9**).
Note: The opposite dummy stays in place.
10. Remove the Dummy Screw.
11. Insert the Clutch Mechanism and torque to 15 Nm (**Fig. 10**).

12. Shorten the aluminum Clutch Mechanism as necessary.
Note: Clutch Key can be used for a more secure grip on the Clutch Mechanism when socket is donned.

Laminated Sockets

1. Put the lower PVA bag over the dummy and tie off with a cord below the upper edge of the dummy. Cut off the excess PVA bag material (**Fig. 1b**).
2. Be careful with the orientation of the Clutch Mechanism. See **figure 5a** for the correct positioning of the black and grey dummies.
Note: The dummies have opposite threads.
3. Use the Central Screw to secure the Lock Body with Silicone Ring on the Manufacturing Dummy and the plaster model. Use the 4-Hole Screws to secure the 4-Hole Dummy with 4-Hole Gasket to the Lock Body. Note: It is important to add the 4-Hole Gasket to prevent resin leaking between Lock Body and 4-Hole Dummy (**Fig. 2b**).
4. Fill the head of the Central Screw and 4-Hole Screws with plasticine or similar. (**Fig. 2b**).
5. Attach layup material according to usual lamination (**Fig. 3b**).
6. Drape unidirectional carbon tape over the distal connector in both A-P and M-L directions and tie off with cord on the base of the Lock Body (**Fig. 4b**).
7. Apply the reinforcement material layup to agree with the expected stress. Tie off with cord on the base of the Lock Body (**Fig. 5b**).
8. Put the upper PVA bag in place and laminate as usual (**Fig. 6b**).
9. After fully hardening, grind and remove the material at the distal end until completely exposing the head of the Central Screw (**Fig. 7**).
10. Remove the Central Screw.
11. Grind carefully at the location of the 4-Hole Dummy until completely exposing the head of the 4-Hole Screws and 4-Hole Dummy (**Fig. 8**).

12. Remove the 4-Hole Screws and 4-Hole Dummy.
13. Grind carefully at the future location of the Clutch Mechanism until completely exposing the head of the Dummy Screw (**Fig. 9**).
Note: The opposite dummy stays in place.
14. Remove the Dummy Screw.
15. Insert the Clutch Mechanism and torque to 15 Nm (**Fig. 10**).
16. Shorten the aluminum Clutch Mechanism as necessary.
Note: Clutch Key can be used for a more secure grip on the Clutch Mechanism when socket is donned.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning. Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 136 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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©Copyright Össur 2021 IFU0643