



Instructions for Use

ICELOCK® 234 LANYARD

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	24
NO	Bruksanvisning	29
DA	Brugsanvisning	34
SV	Bruksanvisning	39
EL	Οδηγίες Χρήσης	44
FI	Käyttöohjeet	49
NL	Gebruiksaanwijzing	54
PT	Instruções de Utilização	59
PL	Instrukcja użytkowania	64
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	74
JA	取扱説明書	79
ZH	中文说明书	84
KO	사용 설명서	89

ENGLISH

MD

Medical Device

DESCRIPTION

The device is a Lanyard Lock, a lock system to be used in conjunction with locking liners.

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.

The weight limit for the device is 120 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Lamination

1. Flatten the Distal End (1) enough to centrally attach the Delrin Guide (2) with the four Nails (3).
2. Apply the inner PVA bag. Pull through the central hole of the Aluminum Ring (4). Flare the bag away from the central hole and attach the ring using the Delrin Screw (5), sealing the PVA bag in the process. Ensure that the Arrow (6) on the Aluminum Ring is aligned in the AP or ML direction. Coat the four Spacing Screws (7) and the Delrin Screw with a release agent.
3. Tie in your lay-up material. Tie all but the first layer beyond the Aluminum Ring. Ensure that the lamination material is at least 6 mm thick (8) at the distal end and the four Spacing Screws are embedded evenly in the material.
4. Apply outer PVA bag and finish lamination.
5. Grind the distal end accurately to expose the Spacer Screw heads.
Warning: Do not grind any deeper.
6. Remove the Spacer Screws and attach the Lock Body (9) with the four long Countersunk Screws (10). Use medium strength threadlocker and torque to 10 Nm.

Thermoplastic Sockets

1. Flatten the Distal End (1) enough to centrally attach the Delrin Guide (2) with the four Nails (3).
2. Starting at the edge of the Delrin Guide, drill a small hole with a wire drill through the full length of your plaster model to its base. This assures maximum vacuum at the distal end.

3. Attach the Aluminum Ring (4) to the Delrin Guide using the Delrin Screw (5). Ensure that the Arrow (6) on the Aluminum Ring is aligned in the AP or ML direction.
4. Vacuum form the thermoplastic. Ensure proper forming around the Aluminum Ring. It is important that the thermoplastic is at least 6 mm thick (8) at the distal end.
5. Grind the distal end accurately to expose the Spacer Screw heads. Do not grind any deeper.
6. Remove the Spacer Screws (7) and attach the Lock Body (9) with the four long Countersunk Screws (10). Use Loctite 410/411 and torque to 10 Nm.

Lanyard Components

1. Shorten the Elbow (11) to the desired length and glue it into the bushing (12) with Loctite 410/411.
2. Find the correct positioning for the Cleat (13), drill two 3 mm holes, and rivet it to the socket.
3. Trim the Cord Cover (14) to the desired length. Glue it into the elbow with Loctite 410/411 and to the socket in the correct location.

Cleat positioning (Fig. C)

(16) Released

(17) Locked

Caution: Tie a knot in the Lanyard Cord to stop it from inadvertently being pulled into the channel.

Caution: It is recommended that the Lanyard Cord (15) be regularly checked for wear.

Note: The Pull Handle can be used for getting a more secure grip on the Lanyard Cord when the socket is donned.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.


DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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