ENGLISH

MD Medical Device

DESCRIPTION

The Icelock 400 for Double Lamination is a mechanical lock that offers a smooth transition and locks in every position. It is easily interchangeable and removable from the inside of the socket.

Components (Fig. 1)

1. Pin
2. Funnel
3. Lock Body
4. Release Mechanism
5. Push Button
6. Protective Sleeve
7. Dust Cap
8. Short Mounting Screws
9. Coupling Plate
10. Long mounting screws
11. Lamination ring

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

INDICATIONS FOR USE AND/TARGET PATIENT POPULATION

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

To fabricate the socket it requires:

- Fixture (P-20109)
- Dummy (P20151)

1. Flatten the distal end of plaster cast and attach the Fixture with two small nails. The transition from cast to the Fixture should be smooth. Ensure the Fixture is mounted in alignment with the residual limb (Fig. 2).
2. Fabricate a soft liner of 3mm PeLite or similar material. Leave the distal bolt end open (Fig. 3).
3. Mount the bolt to the Fixture (Fig. 4).
4. Apply the Dummy on top of the Fixture. Ensure that the bulge for the push button is orientated correctly (Fig. 5).
5. Insulate with PVA-bag and seal at the top of the Dummy (Fig. 6).
6. Place the Lamination Ring on top of the Dummy between the last two layers of reinforcement material. Laminate as usual (Fig. 7).
7. Remove the cured push from the cast. Drill carefully to reveal the center screw in the bulge for the Push Button. Push the plastic piece into the Dummy. Remove Dummy from the socket. Grind a hole for the Push Button (Fig. 8).
8. Attach the Coupling Plate to the Lock Body with the Short Mounting Screws. Then insert the Coupling Plate and Lock Body into the socket (Fig. 9).
9. Drill 3 holes into the Lamination Ring through the Coupling Plate. Secure the Coupling Plate to the Lamination Ring with the Long Mounting Screws (Fig. 10).
10. Trim push button to desired length and screw onto the Lock Body (Fig. 11).
11. Apply the Funnel to the Lock Body from the inside of the socket (Fig. 12).

Attachment Pin

Apply Loctite 410/411 onto threads of Attachment Pin Smooth and screw into the locking liner, torque to 4 Nm with a 7 mm torque wrench.

Caution: Use only Attachment Pin Smooth. The locking pin should be exchanged at least once a year.

Caution:

- Do always check that the threaded adapter for the locking pin is tightened before the prosthesis is attached.
- Maximum weight of the prosthesis: 10 kg.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed. No contact with salt water or chlorinated water is allowed. Dry with a cloth after contact with fresh water or humidity. Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution:

Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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