



Instructions for Use

ICELOCK[®] 551
EXPULSION VALVE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	16
NO	Bruksanvisning	19
DA	Brugsanvisning	22
SV	Bruksanvisning	25
EL	Οδηγίες Χρήσης	28
FI	Käyttöohjeet	31
NL	Gebruiksaanwijzing	34
PT	Instruções de Utilização	37
PL	Instrukcja użytkowania	40
TR	Kullanım Talimatları	43
RU	Инструкция по использованию	46
JA	取扱説明書	49
ZH	中文说明书	52
KO	사용 설명서	55

DESCRIPTION

The device is an auto-expulsion valve with push button for suction release.

Components: (Fig. 1 & 2)

1. Valve Housing
2. Rubber Membrane
3. Cap
4. Cap Key
5. Socket Insert
6. Set Screw
7. Foam Cap

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Directly to the socket

1. Fabricate the socket as usual.
2. Mark the position of the valve onto the socket.
3. Using an 8.5 mm drill bit, make a hole perpendicular to the socket wall.
4. Tap the hole using M10x1.5 tap.
5. Ensure the Cap and Rubber Membrane have been removed prior to shortening the threads to prevent debris lodging in the airways.
6. Measure socket wall thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.

Warning: Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off (**Fig. 3**).

Using the Socket Insert (Fig. 2)

1. The Socket Insert is used as a fixation device for the valve, and should always be used when fabricating a flexible inner socket.
2. Prepare the plaster cast as usual.
3. Fix the Socket Insert to the cast using the Set Screw. Ensure it sits flush to the cast.
4. Plaster may have to be removed or added.
5. Apply the Foam Cap to the Set Screw.
6. Fabricate the socket as usual.
7. Grind material away over the socket insert until Foam Cap is fully exposed.
8. Unscrew the Set Screw to remove socket from cast.
9. Ensure the Cap and Rubber Membrane have been removed prior to shortening the threads to prevent debris lodging in the airways.
10. Measure socket wall and Socket Insert thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.

Warning: Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off (Fig. 3).

Valve Installation

1. Apply silicone adhesive to the threads and housing, where it comes in contact with the socket.
Note: It is not recommended to use epoxy or similar as it can compromise the seal.
2. Install valve and hand tighten with a 5 mm hex key.
Note: Be careful not to over-tighten as it may damage the threads.
Note: After installation distribute silicone adhesive around the valve / socket contact point to ensure an air tight seal.
3. Ensure the valve parts and airways are free from debris before final assembly.
4. Replace Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity. (Fig. 4)
Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.
5. Screw the Cap on using the Cap Key.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Disassembly

1. Unscrew the Cap using the Cap Key.
2. Remove the Rubber Membrane carefully.
3. Clean all parts with water and mild detergent.
4. Hold down the release button while cleaning airways.

Reassembly

1. Replace Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity (**Fig. 4**).
2. Screw the Cap on using the Cap Key.

Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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