## ENGLISH

**MD** Medical Device

- 
- 
-

**DESCRIPTION**

The device is an auto expulsion valve for use in a vacuum socket system.

## Components (Fig. 1)

- Bottom Attachment with Expulsion Valve
- Fabrication Ring with O-Ring, Silicone Ring, and Silicone Seal (Included in L-641000)
- Cross Driver (Included in L-680000)

## INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb. The device must be fitted and adjusted by a healthcare professional. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional. It is intended to be used with a cushion liner and a sleeve.

## Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications. The device is for high impact use, e.g., walking and occasional running.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional. The device is for single patient use.

### ASSEMBLY INSTRUCTIONS

#### Following the socket manufacture

- Make sure that the seal between the Fabrication Ring and Attachment is uncontaminated, without debris.

## LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.

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#### USAGE

**Cleaning and care**
Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

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#### REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

#### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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