

Instructions for Use

ICELOCK® 672 MALE PYRAMID TI
ICELOCK® 673 MALE PYRAMID AL

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	17
ES	Instrucciones para el uso	23
IT	Istruzioni per l'uso	29
NO	Bruksanvisning	35
DA	Brugsanvisning	41
SV	Bruksanvisning	47
EL	Οδηγίες Χρήσης	53
FI	Käyttöohjeet	60
NL	Gebruiksaanwijzing	66
PT	Instruções de Utilização	72
PL	Instrukcja użytkowania	78
TR	Kullanım Talimatları	84
RU	Инструкция по использованию	90
JA	取扱説明書	97
ZH	中文说明书	103
KO	사용 설명서	109

Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

For more detailed device information, see **Table 1**.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

For weight limit and impact level, see **Table 2**.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Components:

- Stainless Steel Pyramid 673 (3)

OR

- Titanium Pyramid 672 (3)

AND

- Pyramid Cap (2)
- Pyramid Bolt w/ O-Ring (1)
- Plastic Washer M8 (1b)
- Silicone Ring (11)
- Silicone Seal (6)

Note: Use with L-680000 Fabrication Tooling Kit (reusable)

When fabrication is complete, insert one of the following locking components.

- L-621000 Icelock ratchet 621.
- L-631000 Icelock lanyard 631.
- L-641000 Icelock valve 641
- L-651000 Icelock smooth 651.

Lamination Instructions:

(not intended for thermoplastic sockets)

1. Flatten distal end of Plaster Cast (9) to center Steel Guide (7). Attach with double-sided tape or two nails maintaining original pin alignment.
2. Apply inner PVA bag and tie it around “neck” of Steel Guide (7).
3. Trim bag around string and place Silicone Seal (6) over “neck” of Steel Guide (7), onto PVA bag.
4. Screw Fabrication Ring (4) (from locking component kit) onto Lock Dummy (5).

5. Choose appropriate length of Release Dummy (8) and screw it firmly into Lock Dummy (5).
6. Apply release agent to components to enable easy removal after lamination (I-400300 Silicone Adhesive is recommended).
7. Place the components onto Steel Guide (7) (make sure that O-Ring on top of Lock Dummy is in place), Pyramid (3) onto Lock Dummy (5), Pyramid Cap (2) and Plastic Washer (1b) onto Pyramid (3) and secure them with Pyramid Bolt (1) (make sure that O-Ring on Bolt is in place).
Note: Be sure to tighten (1-2Nm) so that resin cannot leak in between components.
8. Apply clay into hexagon hole on Release Dummy (B) and Pyramid Bolt (C).
9. Tie lay-up material into grooves on Fabrication Ring (4), Pyramid (3) and Pyramid Cap (2).
Warning: Tie at least two layers of fibers in groove on Pyramid Cap (2).
10. Apply outer PVA bag and pour resin
11. When curing has finished, remove Bolt (1) and Pyramid Cap (2), grind down to Release Dummy (8) and remove.
Warning: Do not grind material away from distal end of socket since that reduces strength of connection between socket and pyramid.
12. Remove socket from plaster cast.
13. Remove Steel Guide (7) and Silicone Seal (6).
14. Use Hex Driver (10) and socket wrench to release Lock Dummy (5).
15. Insert appropriate locking component according to instructions enclosed.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

See which environmental conditions apply to the device in **Table 1**.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Dry with a cloth after contact with fresh water or humidity.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

SPECIFICATIONS

Table 1 - Device Specifications

Part #	Description	Environmental Condition
L-673000	Icelock 673 Stainless Steel Pyramid	Weatherproof
L-672000	Icelock 672 Titanium Pyramid	Waterproof

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

Table 2 - Usage Specifications

Part #	Impact Level	Weight Limit	ISO 10328 - Label Text
L-673000	The device is for low to moderate impact use, e.g., walking.	100 kg	ISO 10328 - P5 - 100 kg
L-672000	The device is for low to high impact use, e.g., walking and occasional running.	166 kg	ISO 10328 - P7 - 166 kg

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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©Copyright Össur 2022-05-10