


Instructions for Use

ICELOCK® 214 4-HOLE CLUTCH

EN		Instructions for Use	4
DE		Gebrauchsanweisung	9
FR		Notice d'utilisation	15
ES		Instrucciones para el uso	21
IT		Istruzioni per l'uso	27
NO		Bruksanvisning	33
DA		Brugsanvisning	39
SV		Bruksanvisning	45
EL		Οδηγίες Χρήσης	51
FI		Käyttöohjeet	57
NL		Gebruiksaanwijzing	63
PT		Instruções de Utilização	69
PL		Instrukcja użytkowania	75
CS		Návod k použití	81
TR		Kullanım Talimatları	87
RU		Инструкция по использованию	93
JA		取扱説明書	99
ZH		中文说明书	104
KO		사용 설명서	109

1

ENGLISH

Medical Device

DESCRIPTION

The device is a mechanical lock used to connect and release a locking liner into a prosthetic socket.

Components (Fig. 1)

1. Delrin Screw
2. Delrin Tooling Body
3. Plaster Cast
4. Lam/Therm Ring
5. Lock Body
6. Countersunk Screws
7. Clutch Mechanism
8. Push Button
9. T-Wrench
10. Attachment Pin

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.

The weight limit for the device is 120 kg.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Lamination Instructions

1. Flatten distal end of Plaster Cast enough to centrally attach Delrin Tooling Body. Attach with nails provided or double-sided tape maintaining original pin alignment.
2. Apply inner PVA bag. Pull through central hole of Lam/Therm Ring. Flare bag away from the central hole and attach ring using Delrin Screw, sealing PVA bag in the process. **Note:** Ensure that arrow (A) on Lam/Therm Ring is aligned in AP or ML direction.
3. Place some clay in spacing screw heads (B) and in a groove on Delrin Screw (C).
4. Tie in lay-up material. Tie all but first layer beyond Lam/Therm Ring. Ensure that four spacing (B) screws are embedded evenly in material.
5. Apply outer PVA bag and complete lamination.

6. Grind distal end accurately to expose spacer screw heads (B).
Warning: Do not grind any deeper.
7. Remove spacer screws (B) and attach Lock Body with four long Countersunk Screws. Use medium strength threadlocker and torque to 10 Nm.

Note: Push Button may be shortened by unscrewing from Lock Body and cutting it with hacksaw.

Thermoplastic Instructions

1. Flatten distal end of Plaster Cast enough to centrally attach Delrin Tooling Body. Attach with nails provided or double-sided tape maintaining original pin alignment.
2. A small hole may be drilled from edge of Delrin Tooling Body (D) through length of Plaster Cast to assure maximum vacuum.
3. Attach Lam/Therm Ring to Delrin Tooling Body with Delrin Screw.
Note: Ensure that arrow (A) on Lam/Therm Ring is aligned in AP or ML direction.
4. Vacuum form thermoplastic. Ensure proper forming, especially at distal end. Thermoplastic must be at least 6 mm thick at distal end.
5. Rough-cut and grind plastic in the transverse plane, stopping at spacer screws (B). **Warning:** Do not grind any deeper.
6. Remove spacer screws (B) and attach Lock Body with four long Countersunk Screws. Use medium strength threadlocker and torque to 10 Nm.

Note: Push Button may be shortened by cutting it with hacksaw.

Warning: If Clutch Mechanism is removed during lifetime of lock, torque to 15 Nm.

Attachment Pin

Apply Loctite 410/411 onto the threads on the Attachment Pin and screw it into the Locking Liner. Torque the Attachment Pin to 4 Nm.

Caution: Use Attachment Pins for Icelock 200 series only.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning. Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.

Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P5 - 120 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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Össur UK Ltd