


Instructions for Use

ICELOCK® 674 4-HOLE ADAPTER AL
ICELOCK® 675 4-HOLE ADAPTER TI

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	17
ES	Instrucciones para el uso	24
IT	Istruzioni per l'uso	31
NO	Bruksanvisning	38
DA	Brugsanvisning	44
SV	Bruksanvisning	50
EL	Οδηγίες Χρήσης	56
FI	Käyttöohjeet	63
NL	Gebruiksaanwijzing	69
PT	Instruções de Utilização	76
PL	Instrukcja użytkowania	83
TR	Kullanım Talimatları	90
RU	Инструкция по использованию	96
JA	取扱説明書	103
ZH	中文说明书	109
KO	사용 설명서	114

ENGLISH

MD

Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

For more detailed device information, see **Table 1**.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

For weight limit and impact level, see **Table 2**.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur replacement parts for the device for ensured durability

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Components:

- Aluminum 4-Hole 674 (3)
- Screws M6x20
- Silicone Seal (5)

OR

- Titanium 4-Hole 675 (3)
- Screws M6x20 Heavy Duty
- Silicone Seal (5)

Note: Use with L-680000 Fabrication Tooling Kit (reusable)

Note: When fabrication is complete, insert one of the following locking components (include Fabrication Ring (4) needed in the fabrication process)

- L-621000 ICELOCK ratchet 621
- L-631000 ICELOCK lanyard 631
- L-641000 ICELOCK valve 641
- L-651000 ICELOCK smooth 651

Lamination Instructions:

1. Flatten distal end of Plaster Cast (8) to center Steel Guide (6).
2. Attach with double-sided tape or two nails maintaining original pin alignment.
3. Apply inner PVA bag and tie it around “neck” of Steel Guide (6).
4. Trim bag around string and place Silicone Seal (5) over “neck” of Steel Guide (6), onto PVA bag.
5. Screw 4-Hole Adapter (3) onto Lock Dummy (4).
6. Choose appropriate length of Release Dummy (7) and screw it firmly into Lock Dummy (4).

7. Screw Pin Dummy Bolt (2) into shorter Pin Dummy (1) making sure that hexagon hole on Pin Dummy Bolt stays visible.
8. Place the components onto Steel Guide (6) and secure them with Pin Dummy and Bolt (1;2) (make sure that O-Ring on top of Lock Dummy is in place).
Note: Be sure to tighten (1-2 Nm) so that resin cannot leak in between components.
9. Apply release agent to components (A) to enable easy removal after lamination (I-400300 Silicone Adhesive is recommended).
10. Apply clay into hexagon hole on Release Dummy (B) and through hole on Pin Dummy (C).
11. Tie lay-up material into groove on 4-Hole Adapter (3) and over components.
Warning: Tie at least two layers of fibers in groove and over the distal part of adapter.
12. Apply outer PVA bag and pour resin.
13. When curing has finished, grind down to Release Dummy (7) and remove.
14. Remove socket from plaster cast.
15. Remove Bolt (2) in Pin Dummy (1) from inside of socket.
16. Remove Steel Guide (6) and Silicone Seal (5).
17. Use Hex Driver (9) and socket wrench to release Lock Dummy (4).
18. Remove Pin Dummy (1).
19. Grind down to spacers on 4-Hole Adapter (3) and remove all four silicone sealings.
Warning: Do not grind any further since that reduces strength of connection between socket and adapter.
20. Insert appropriate locking component according to instructions enclosed.

Warning: For secure fastening use medium strength threadlocker and apply rated torque twice, see **Table 1**.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

See which environmental conditions apply to the device in **Table 1**.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

A Waterproof device can be used in a wet or humid environment and submerged in up to 1-meter-deep fresh water for a maximum of 30 minutes. No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

SPECIFICATIONS

Table 1 - Device Specifications

Part #	Description	Tightening Torque	Environmental Condition
		Mont Screws	
L-674000	Icelock 674 Aluminum 4-Hole	10 Nm	Weatherproof
L-675000	Icelock 675 Titanium 4-Hole		Waterproof

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

Tale 2 - Usage Specifications

Part #	Impact Level	Weight Limit	ISO 10328 - Label Text
L-674000	The device is for low to moderate impact use, e.g., walking.	100 kg	ISO 10328 - P5 - 100 kg
L-675000	The device is for low to high impact use, e.g., walking and occasional running.	166 kg	ISO 10328 - P7 - 166 kg

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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Össur UK Ltd