THE ÖSSUR NEWS

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Life Without Limitations

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1. Kristinsson Ö. The ICEROSS concept: a discussion of a philosophy. *Prosthet Orthot Int* 1993;17:49-55.
2. Narita H, Yokogushi K, Skii S, Kakizawa M, Nosaka T. Suspension effect and dynamic evaluation of the total surface bearing (TSB) trans-tibial prosthesis: a comparison with the patellar tendon bearing (PTB) trans-tibial prosthesis. *Prosthet Orthot Int* 1997;21:175-178.
3. Karp J. New technologies, materials improve patient options. *BioMechanics Magazine Online*.
4. Kahle, JT. Conventional and hydrostatic transtibial interface comparison. *J. Prosthet Orthot* 1999;11(4):85-91.
5. Cluitmans J, Geboers M, Deckers J, Rings F. Experiences with respect to the Iceross system for trans-tibial prostheses. *Prosthet Orthot Int* 1994;18:78-83.
6. Larsson GV, Johannesson A, Holmquist A, Larsson B. Lower extremity amputations: a controlled active treatment protocol. *Orthopediskt Magasin* 1997;4:20-21.
7. Datta D, Vaidya SK, Howitt J, Gopalan L. Outcome of fitting an Iceross prosthesis: views of trans-tibial amputees. *Prosthet Orthot Int* 1996;20:111-115. Yigiter K, Sener G, Bayar K. Comparison of the effects of patellar tendon bearing and total surface bearing sockets on prosthetic fitting and rehabilitation. *Prosthet Orthot Int* 2002;26:206-212.
8. Tazawa E. The TSB socket and silicone suction sleeve concept. *Jap J Prosthet Orthot* 1991; 7:145-150.
9. Dasgupta AK, McCluskie PJA, Patel VS, Robins L. The performance of the ICEROSS prosthesis amongst transtibial amputees with a special reference to the workplace- a preliminary study. *Occup. Med.* 1997; 47:228-236.
10. Winn, Alastair. Factors in selecting medical silicones. *Medical Plastics and Biomaterials* 1996;3(2): 16-19.
11. Heide C. Silicone rubber for medical device application. *Medical Device and Diagnostic Industry* November 1999 column.
12. Heim M, Wershavski M, Zwas ST, Siev-Ner I, Nadvorna H, Azaria M. Silicone suspension of external prostheses. *J Bone Joint Surg [Br]* 1997;79-B:638-40.
13. Fillauer CE, Pritham CH, Fillauer KD. Evaluation and development of the silicone suction socket (3S) for below-knee prostheses. *J Prosthet Orthot Int* 1989;1:92-103.
14. Postema K, Hermens HJ, de Vries J, Koopman HF, Eisma WH. Energy storage and release prosthetic feet. Part 2: Subjective ratings of 2 energy storing and 2 conventional feet, user choice of foot and deciding factor. *Prosthet Orthot Int* 1997;21(1):28-34.
15. Cortes A, Viosca E, Hoyos JV, Prat J, Sanchez-Lacuesta J. Optimization of the prescription for trans tibial (TT) amputees. *Prosthet Orthot Int* 1997;21(3):168-74.
16. Day HBJ. The Knud Jansen Lecture 1998. Amputation rehabilitation – finding the niche.
17. Davies B, Datta D. Mobility outcome following unilateral lower limb amputation. *Prosthet Orthot Int* 2003;27:186-190.
18. Hoaglund FT, Jergesen HE, Wilson L, Lamoreux LW, Roberts R. Evaluation of problems and needs of veteran lower-limb amputees in the San Francisco Bay area during the period 1977-80. *J Rehabil Res Dev* 1983;20(1):57-71.
19. Kegel B, Webster JC, Burgess EM. Recreational activities of lower extremity amputees: a survey. *Arch Phys Med Rehabil* 2000;81(8):1039-44.
20. Legro MW, Reiber G, del Aguila M, Ajax M, Boone DA, Larsen JA, Smith DG, Sangeorzan B. Issues of importance reported by persons with lower limb amputation and prostheses. *J. Rehabil Res Dev* 1999;36(3):155-63.
21. Nielsens CC. A survey of amputees: functional level and life satisfaction, information needs, and the prosthetist's role. *J Prosthet Orthot* 1990;3:125-9.
22. Pell JP, Donnan PT, Fowkes FG, Ruckley CV. Quality of life following lower limb amputation for peripheral arterial disease. *Eur J Vasc Surg* 1993;7:448-451.
23. Pitkin MR. Effects of design variants in lower-limb prostheses on gait synergy. *J Prosthet Orthot* 1997;9(3):113-122.
24. Pohjolainen T, Alaranta H, Karkainen M. Prosthetic use and functional and social outcome following major lower limb amputation. *Prosthet Orthot Int* 1990;14:75-79.
25. Klute GK, Kallfelz CF, Czerniecki JM. Mechanical properties of prosthetic limbs: adapting to the patient. *J Rehabil Res Dev* 2001;38(3):299-307.
26. Falstie-Jensen N, Christensen K-S, Brochner-Mortensen J. Long posterior flap versus equal sagittal flaps in below-knee amputation. *J-Bone-Joint-Surg-[Br]* 1989;71(1):S102-4.
27. Gailey R. Stability within the socket creates a stable world. *O&P Edge*. September 2003.
28. Murdoch G. The Knud Jansen Lecture 1983. Amputation Revisted.
29. Murdoch G. Myoplastic Techniques. *Bull Prosthet Res* 1968;10-9:4-13.
30. Ruckley CV, Stonebridge PA, Prescott RJ. Skewflap versus long posterior flap in below-knee amputations: multicenter trial. *J Vasc Surg* 1991;13(3):423-7.
31. Stokosa JJ. New Developments in Prosthetics, In: Moore WS and Malone JM, eds. *Lower Extremity Amputation*, W.B. Saunders Company, 1989.
32. Stokosa JJ. Prosthetics for Lower Limb Amputees, In: Haimovici H, Ascer E, et al eds. *Haimovici's Vascular Surgery*, 4th ed, Boston: Blackwell Science, Inc, 1996.
33. Walter Reed Medical - Center for Amputation Surgery and Rehabilitation with the Amputee Coalition of America Conference on Amputation Surgical Techniques – March 11, 2003.
34. Zhang M, Fan Y. Biomechanical studies on lower-limb prosthetic socket design. *Advances in biomechanics*, 2001, Beijing, China 134-38.
35. Murdoch G. The Dundee socket – a total contact socket for the below-knee amputation. *Health Bulletin*, Vo. XXII: No. 4, pp. 70-71. October 1968.
36. Radcliffe CW, Foort J. The patellar-tendon-bearing below-knee prosthesis. *Biomechanics Laboratory*, University of California, Berkeley and San Francisco, 1961.
37. Wilson AB. Evaluation of the patellar-tendon-bearing prosthesis and its variations. *Prosthetic and Orthotic Practice* (ed. George Murdoch). E. Arnold, London, 1969.

Life Without Limitations

38. Sewell P, Noroozi S, Vinney J, Andrews S. Developments in the trans-tibial prosthetic fitting process: a review of past and present reseach. *Prosthet Orthot Int* 2000; 24:97-107.
39. Levy S. Skin probles of the leg amputee. *Prosthet Orthot Int* 1980;4:37-44. Goldstein B, Sanders J. Skin response to repetitive mechanical stress: a new experimental model in pig. *Arch Phys Med Rehabil* 1998;79(3):265-72.
40. Lyon CC, Kulkarni J, Zimerson E, Van Ross E, Beck MH: Skin disorders in amputees. *J Am Acad Dermatol* 2000;42:501-7.
41. Sonck WA, Cockrell JL, Koepke GH. Effect of liner materials on interface pressures in below-knee prostheses. *Arch Phys Med Rehabil Nov.* 1970:666-669.
42. Zhang M, Turner-Smith AR, Tanner A, Roberts VC. Frictional action at residual limb/prosthetic socket interface. *Med Eng & Phys* 1996;18(3):207-14.
43. Zhang M, Turner-Smith AR, Roberts VC. The reaction of skin and soft tissue to shear forces applied externally to the skin surface. *J of Eng in Medicine* 1994;208(4):217-22.
44. Buis AWP, Concery P. Socket/stump interface dynamic pressure distributions recorded during the prosthetic stance phase of gait of a trans-tibial amputee wearing a hydrocast socket. *Prosthet and Orthot Int* 1999;23:107-12.
45. Convery P, Buis AWP. Conventional patellar-tendon-bearing socket stump interface dynamic pressure distributions recorded during prosthetic stance phase of gait of a trans-tibial amputee. *Prosthet Orthot Int* 1998;22(3):193-98.
46. Sanders JE, Daly CH. Interface pressures and shear stresses: sagittal plane angualr alignment effects in three trans-tibial amputee case studies. *Prosthet Orthot Int* 1999;23:21-29.
47. Sanders JE, Lam D, Dralle AJ, Okumura R. Interface pressures and shear stresses at thirteen socket sites on two persons with transtibial amputation. *J. Rehabil Res Dev* 1997;34(1):19-43.
48. Sanders JE, Daly CH, Burgess EM. Interface shear stresses during ambulation with a below-knee prosthetic limb. *J Rehabil Res Dev* 1992;29(4)1-8.
49. Zhang M, Turner-Smith AR, Tanner A, Roberts VC. Clinical investigation of the pressure and shear stress on the trans-tibial residuum with prosthesis. *Med Eng Phys* 1998;20:188-98.
50. Erikson U, Lempberg R. Roentgenological study of movements of the amputation stump within the prosthesis socket in below knee amputees fitted with a PTB prosthesis. *Acta Orthopaedica Scandinavica* 1969;40:520-9.
51. Grevsten S, Erikson U. A roentgenological study of the stump-socket contact and skeletal displacement in the PTB-suction prosthesis. *Uppsala Journal of Medicine Science* 1975;80:49-57.
52. Sabolich J, Guth T. Below-knee prosthesis with total flexible socket (t.f.s): a preliminary report. *Clinic Prosthet and Orthot* 1986;10(2):93-99.
53. Tanner JE, Berke GM. Radiographic comparison of vertical tibial translation using two types of suspensions on a transtibial prosthesis: a case study. *J Prosthet Orthot* 2001;13(1):14-16.
54. Wirta RW, Golbranson FL, Mason R, Calvo K: Analysis of below-knee suspension systems: effects on gait. *J Rehabil Res Dev* 1990;27:385-396.
55. Narita H, Yokogushi K, Skii S, Kakizawa M, Nosaka T. Suspension effect and dynamic evaluation of the total surface bearing (TSB) trans-tibial prosthesis: a comparison with the patellar tendon bearing (PTB) trans-tibial prosthesis. *Prosthet Orthot Int.* 1997;21:175-178.
56. Lake C, Supan T. The incidence of dermatological problems in the silicone suspension sleeve. *J Prosthet Orthot* 1997;9(3):97-106.
57. Smith MW. Effects of aging on skin. WebMD Medical Reference in collaboration with The Cleveland Clinic. April 2003.
58. The skin from Seagull Cosmetics website. (<http://www.seagullcosmetics.com>)
59. Zhang YP, Mak AFT. In vivo properties of human skin. *Prosthet Orthot Int* 1999;23:135-41.
60. Boughton B. Comorbidities complicate management of amputees. *BioMechanics Magazine Online.*
61. Stephens E. Peripheral Vacular Disease. WebMD Medical Reference in collaboration with The Cleveland Clinic. April 2003.
62. uabhealthsystem website: www.health.uab.edu
63. Hampton Roads Health System website: www.bonsecourshamptonroads.greystone.net
64. International Diabetes Federation website: www.idf.org
65. Consensus development conference on diabetic foot wound care. *Diabetes Care* 1999; 22: 1354–60.
66. Diabetes Atlas, second edition, International Diabetes Federation, 2003.
67. Diabetes and Cardiovascular Disease: Time to Act, International Diabetes Federation, 2001.
68. Website of the World Health Organisation. Diabetes around the world, International Diabetes Federation, 1998.
69. Cost-effective Approaches to Diabetes Care and Prevention, International Diabetes Federation, 2003.
70. Board WJ, Street GM, Caspers C. A comparison of transtibial amputation suction and vacuum assisted socket conditions. *Prosthet Orthot Int* 2001;25:202-9.
71. Chino N, Pearson JR, Cockrell JL. Negative pressure during swing phase in below-knee prostheses with rubber sleeve suspension. *Arch Phys Med Rehabil* 1975;56:22-26.
72. Grevsten S. Patellar Tendon Bearing Suction Prosthesis. *Upsala J. Med Sci* 1977;82:209-220.
73. Grevsten S, Marsh L. Suction type prosthesis for below knee amputees, a preliminary report. *Artif Limbs* 1971;15:78-80.
74. Infield S. Key benefits of suction socket suspension systems for below-kjnee prostheses. *O&P Business News*, June 1, 2001 pp. 32-34.
75. Murphy E. Lower extremity components. *Orthopaedic Appliance Atlas*, Vol 2, JW Edwards, 1960.
76. Pearson JR, Grevsten S, Almby B, Marsh L. Pressure variation in the below-knee patellar tendon bearing suction socket prosthesis. *J Biomechanics* 1974;7:487-96.
77. Staats TB, Lundt J. The UCLA total surface bearing suction below-knee prosthesis. *Clin Prostet Orthot* 1987;11(3):118-30.
78. Össur R&D. Internal paper regarding transport of active (skin care) materials in second generation liners.

Life Without Limitations